Rose Water Facial Toner

with Hyaluronic Acid



You can maintain these proper functions with the

Rose Water Facial Toner:







Packaging: Clear bottle with spray dispenser

Size: 4.0 fl oz

Dimensions: 2" x 2" x 6.5" Gross Weight: 3.6 fl oz







Free



Free



Rose Water Facial Toner

Soothing and Calming Toner

TARGET MARKET:

Consumers that are looking for a facial toner that moisturizes and soothes the skin.



Directions:

Spritz face and neck area after cleansing. Apply as needed throughout the day. May be applied on top of makeup. For best results, follow with your favorite moisturizer.

Warnings:

Store at room temperature. For external use only. Keep out of reach of children.

Ingredients:

Rosa Damascena Rose Flower Water (Hydrosol), Aloe Barbadensis Leaf Juice, Vegetable Glycerin, Hamamelis Virginiana (Witch Hazel), Water, Phenoxyethanol, Hyaluronic Acid (Sodium Hyaluronate).

Hyalogic's Rose Water Facial Toner is made with steamed distilled rose water and Hyaluronic Acid to hydrate, soothe and leave your skin feeling fresh and awake. This toner is great for prepping your skin for a moisturizer, as an excellent makeup refresher or as a midday pick-me-up.

Hyaluronic Acid (HA) is a natural and vital element of healthy, supple skin. HA has the unique ability to bind up to 1,000 times its weight in water, which makes it nature's greatest natural hydrator for your skin!

Rose Water is created by distilling rose petals with steam. Rose water is fragrant, and it's sometimes used as a mild natural fragrance as an alternative to chemical-filled perfumes. Rose water has been used for thousands of years, including in the Middle Ages. It's been used traditionally in both beauty products and food and drink products. It also comes with plenty of potential health benefits, including high in Vitamins A and C which are known to rejuvenate the skin.



Vancouver BC 1.800.665.7065 orders@ecotrend.ca www.ecotrend.ca



Manufactured for:

610 NW Platte Valley Dr. Riverside, MO 64150 USA Toll-Free: 866-318-8484 Website: www.hyalogic.com **REV 072718**